
Anti-PD-L1 / DIA-PDL1-OD

Mouse monoclonal anti-cancer cell marker

Clone JAL1

Product Information

Catalog No.:	DIA-PDL1-OD	Presentation:	Purified antibody in Tris pH 7.3-7.7 with 1% BSA, <0.1% NaN ₃
Clone:	JAL1	Applications:	Immunohistochemistry (IHC), standard formalin-fixed paraffin sections
Isotype	Rabbit IgG	Dilutions:	1:100 - 1:200 IHC-P
Quantity	100µl		(General recommendation, validation of antibody performance/protocol is the responsibility of the end user. Positive/negative controls should be run simultaneously with patient specimen. Interpretation must be made by a qualified pathologist within the context of patient's clinical history/other diagnostic tests.)
Specificity:	PD-L1		
Physical State:	Liquid		
Species			
Reactivity:	Human		
Positive Control:	Tonsil		
Visualization:	Membranous		

Reactivity

Clone JAL1 has been developed specifically for routine immunohistochemical (IHC) detection of PD-L1 in formalin-fixed paraffin-embedded tissue specimen. Moreover, JAL1 has been validated for the identification of PD-L1-positive macrophages and tumor tissues under pathological conditions.

PD-L1, also known as CD274 or B7 homolog 1 (B7-H1) is highly expressed in the heart, skeletal muscle, placenta and lung and weakly expressed in the thymus, spleen, kidney and liver. PD-L1 is expressed on macrophages and activated T- and B-cells, dendritic cells, keratinocytes and monocytes.

Moreover, several human cancer cells express PD-L1 at high levels. Binding of PD-L1 with its receptor PD-1 on T cells inhibits T cell proliferation and the production of cytokines, such as interleukin (IL)-2. It has been shown that PD-L1 helps tumor cells to evade anti-tumor immunity and blockade of PD-L1 reduces the growth of tumors in the presence of immune cells.

PD-L1 is commonly over expressed on tumor cells or on non-transformed cells in the tumor microenvironment². PD-L1 expressed on the tumor cells binds to PD-1 receptors on the activated T cells, which leads to the inhibition of the cytotoxic T cells. These deactivated T cells remain inhibited in the tumor microenvironment.

PD-1 and PD-L1/PD-L2 belong to the family of immune checkpoint proteins that act as co-inhibitory factors, which can halt or limit the development of the T cell response.

Instructions for Use

Immunohistochemical staining of standard formalin-fixed paraffin sections

Deparaffinize and rehydrate according to standard procedures. Heat induced epitope retrieval (HIER) is required (pH 9-10 for 10-30 minutes). For immunohistochemical detection different techniques can be used: indirect immunoenzyme labeling with a secondary antibody conjugate, biotin/(strept)avidin-based detection, soluble enzyme immune complex or polymer-based detection. The antibody can be adapted for use on automated staining instruments.

Intended use / regulatory status

Europe: For in Vitro Diagnostic Use / All other countries: For Research Use only

Storage and Stability

Store at 2-8°C. Do not freeze. The antibody is stable until the date indicated on the label, when stored properly.

Safety Notes


The material contains <1% sodium azide as preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material. Avoid skin and eye contact, inhalation and ingestion.



Figures

Immunohistochemistry of human PD-L1 in routine formalin-fixed paraffin-embedded tissue samples

- A:** Strong PD-L1 positivity in >90% of the cells of a thyroid adenoma.
B: Strong focal PD-L1 positivity in a kidney oncocytoma (mosaic pattern).
C: Thymoma with strong PD-L1 expression in epithelial cells
D: PD-L1 positive Hodgkin's lymphoma with predominant staining of Hodgkin and Reed-Sternberg cells.



(pictures courtesy of Prof. Guido Sauter, Department of Pathology, University Hospital Eppendorf, Hamburg, Germany)

References

1. O'Malley DP et al. Immunohistochemical detection of PD-L1 among diverse human neoplasms in a reference laboratory: observations based upon 62,896 cases. *Modern Pathology* (2019) 32:929–942
2. Büttner R et al. Programmed death-ligand 1 immunohistochemistry testing: a review of analytical assays and clinical implementation in non-small-cell lung cancer. *J Clin Oncol.* (2017) 35:3867–3876
3. Rimm DL et al. A prospective, multi-institutional, pathologist-based assessment of 4 immunohistochemistry assays for PD-L1 expression in non-small cell lung cancer. *JAMA Oncol.* (2017) 3:1051–1058
4. Milne CP et al. Complementary versus companion diagnostics: apples and oranges? *Biomark Med.* (2015) 9:25–34
5. Ott PA et al. CTLA-4 and PD-1/PD-L1 blockade: new immunotherapeutic modalities with durable clinical benefit in melanoma patients. *Clin Cancer Res.* (2013) 19:5300–5309
6. Pardoll DM. The blockade of immune checkpoints in cancer immunotherapy. *Nat Rev Cancer.* (2012) 12:252–264.
7. Yokosuka T et al. Programmed cell death 1 forms negative costimulatory microclusters that directly inhibit T cell receptor signaling by recruiting phosphatase SHP2. *J Exp Med.* (2012) 209:1201–1217

